



"Please support lasing for a cause."

A PARADIGM SHIFT IN THE TREATMENT OF SOFT TISSUE

Gone are the days of using a scalpel or electrosurge. Lasers are the modality of choice for cutting and coagulating tissue, offering considerably less trauma, less bleeding, and less necrosis of tissue. No need to waste time packing cord or worrying about tissue shrinkage. Lasers are the safer, easier way to uncover implants, brackets, and unerupted teeth. Save time, and save money.

"Last year when I introduced Picasso? the world's first dental laser under \$5,000? it was clear that dentists had been waiting for an affordable option. Picasso instantly became the #1 selling laser in the world! Now with the introduction of Picasso Lite at only \$2,495, every practice can put a laser in each operatory? just like other handpieces. Picasso Lite is not only the most affordable laser in the world? it's also the easiest to operate!"


"Welcome to better dentistry and better patient comfort!"

Alan Miller, President/CEO AMD LASERS

THE WORLD'S EASIEST LASER TO OPERATE:

- STEP 1 – Insert disposable tip into the handpiece
- STEP 2 – Select SPEED or COMFORT or DECONTAMINATION
- STEP 3 – Select ready (the red button) and step on the foot control to activate

| | | |
|------------------------------|---|--|
| SPEED: | Surgery Troughing Frenectomies | |
| SUPER PULSE COMFORT MODE: | Gingivectomies Exposing implants & brackets/teeth | |
| DECONTAMINATION: | Pockets Aphthous ulcers Pulpotomies | |



- #1 Selling Laser in the World.
- Top Product 2009
- Endorsed and Loved by Top Lecturers, Luminaries, Universities, Industry Experts, and Patients Worldwide
- Available in 50+ Countries and Counting
- #1 Laser Training and Support Company in the World
- World's Easiest Laser to Operate
- World's First Laser with Disposable Tips or a Strippable Fiber
- Laser Training and Certification Included

NEW **Picasso** BY AMD LASERS™, LLC
ST **LITE**
\$2,495



- FEATURES:
- 2.5 watts of power
 - Continuous / pulse modes
 - 3 customizable presets
 - Adjustable aiming beam
 - 2-year warranty



A Limited Edition (LE) Picasso Lite is available in pink.

AMD LASERS proudly announces that we are "lasing for a cause" in the fight against breast cancer. Through June of 2010, we have pledged to donate \$15,000 to the National Breast Cancer Foundation.



Picasso BY AMD LASERS™, LLC
\$4,995



- FEATURES:
- 7 watts of power
 - Continuous / pulse modes
 - 8 customizable presets
 - Multiple language menu
 - 3-year warranty



AMD LASERS™, LLC
www.amdlasers.com
7405 Westfield Blvd.
Indianapolis, IN 46240, USA

866.WWW.AMDL
866.999.2635
TEL: +1 317.202.9530
FAX: +1 678.868.4108
sales@amdlasers.com



AMD LASERS™
Effective. Practical. Affordable. World Class Lasers.

ISO 13485